



International Journal of Chemistry and Pharmaceutical Sciences

Journal Home Page: www.pharmaresearchlibrary.com/ijcps



RESEAECH ARTICLE

Method Development and Validation of Elbasvir and Grazoprevir by using RP-HPLC

Sri Mounika, M. Shymala*

Department of Pharmaceutical Analysis, Joginapally B.R. Pharmacy College, Yenkapally, Hyderabad, Telangana-500075

ABSTRACT

A new method was established for simultaneous estimation of Elbasvir and Grazoprevir by RP-HPLC method. The chromatographic conditions were successfully developed for the separation of Elbasvir and Grazoprevir by using Symmetry (4.6 x 250mm, 5 μ m), flow rate was 1.0ml/min, detection wave length was 264nm. The instrument used was WATERS HPLC Auto Sampler, Separation module 2695, PDA Detector, Empower-software version-2. The retention times were found to be 2.8 mins and 4.3 mins. The assay of Elbasvir and Grazoprevir was performed with tablets and the % assay was found to be 100.18 and 100.00 which shows that the method is useful for routine analysis. The linearity of Elbasvir and Grazoprevir was found to be linear with a correlation coefficient of 0.999 and 0.999, which shows that the method is capable of producing good sensitivity. The acceptance criteria of precision is RSD should be not more than 2.0% and the method show precision 0.5 and 0.1 for Elbasvir and Grazoprevir which shows that the method is precise. The acceptance criteria of intermediate precision is RSD should be not more than 2.0% and the method show precision 0.6 and 0.2 for Elbasvir and Grazoprevir which shows that the method is repeatable when performed in different days also. The accuracy limit is the percentage recovery should be in the range of 97.0% - 103.0%. The total recovery was found to be 100.11% and 100.38% for Elbasvir and Grazoprevir. The validation of developed method shows that the accuracy is well within the limit, which shows that the method is capable of showing good accuracy and reproducibility. The acceptance criteria for LOD and LOQ is 3 and 10. The LOD and LOQ for Elbasvir was found to be 2.97 and 9.91 and LOD and LOQ for Grazoprevir was found to be 3.09 and 10.07. The robustness limit for mobile phase variation and flow rate variation are well within the limit, which shows that the method is having good system suitability and precision under given set of conditions.

Keywords: Elbasvir, Grazoprevir, HPLC

ARTICLE INFO

CONTENTS

1. Introduction.	10
2. Materials and Method.	10
3. Results and Discussion.	10
4. Conclusion.	15
5. References.	15

Article History: Received 29 November 2018, Accepted 28 December 2017, Available Online 27 January 2018

*Corresponding Author

M. Shymala
Department of Pharmaceutical Analysis,
Joginapally B.R. Pharmacy College,
Yenkapally, Hyderabad, Telangana-500075
Manuscript ID: IJCPS3514



PAPER-QR CODE

Citation: M. Shymala. Method Development and Validation of Elbasvir and Grazoprevir by using RP-HPLC. *Int. J. Chem, Pharm, Sci.*, 2018, 6(1): 09-15.

Copyright© 2018 M. Shymala. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution and reproduction in any medium, provided the original work is properly cited.

1. Introduction

Analytical methods: Methods are developed for new products when no official methods are available. Alternate methods for existing (non-pharmacopoeial) products are developed to reduce the cost and time for better precision and ruggedness [1,2]. Trial runs are conducted, method is optimized and validated. When alternate method proposed is intended to replace the existing procedure comparative laboratory data including merit/demerits are made available

Description of the Various Analytical Methods

Titrimetric and gravimetric method of analysis is suitable when the sample is present in pure form or when no interference is observed in the mixture with other materials [3]. Ultraviolet and visible spectrometric method is suitable when no Interference is observed in the mixture [4]. HPLC and GC methods are more advantageous than the above due to their capability in separating organic mixtures and quantitative estimations. AAS is used mainly for quantitative estimation in ppm and ppb levels of elements [5]. Infra-red spectroscopy though mainly used for qualitative analysis can be used for quantitative estimation also. Out of all the above methods, thin layer chromatography plays a very important role in analysis due to its adaptability, flexibility, and cost and time. It can be used both for qualitative and quantitative determination. After separation spots can be scanned with the help of a scanner and quantitative measurement can be made [6].




Figure 1: Elbasvir




Figure 2: Grazoprevir

2. Materials and Methods

Apparatus: The instrument used for the study was WATERS, software: Empower, 2695 separation module, UV detector [10].

International Journal of Chemistry and Pharmaceutical Sciences

Preparation of mobile phase:

Accurately measured 450 ml (45%) of above buffer and 550 ml of Acetonitrile HPLC (55%) were mixed and degassed in an ultrasonic water bath for 10 minutes and then filtered through 0.45 μ filter under vacuum filtration.

Preparation of Phosphate buffer:

3.4g of KH_2PO_4 in 1000 ml of HPLC water Ph was adjusted with OPA up to 3.0. final solution was filtered through 0.44 μ m Membrane filter and sonicate it for 10 mins[11].

Diluent: The Mobile phase was used as the diluent.

Optimization Chromatographic trials for Simultaneous Estimation of Elbasvir and Grazoprevir by RP- HPLC.

Optimization chromatographic conditions

Instrument used : Waters HPLC with auto sampler and UV detector.

Temperature : Ambient

Column : Symmetry (4.6 x 250mm, 5 μ m)

Buffer : 3.4g of KH_2PO_4 in 1000 ml of HPLC water Ph was adjusted with OPA up to 3.0.

pH : 3.0

Mobile phase : 45% buffer 55% Acetonitrile

Flow rate : 1 ml per min

Wavelength : 264 nm

Injection volume : 20 μ l

Run time : 10 min.




Figure 3: Optimization Chromatogram

Observation: The separation of two analytical peaks was good. The plate count also above 2000, tailing factor below 2, and the resolution is above 2. The condition is taken as optimized method.

3. Results and Discussion

Method Validation Parameters

1. Specificity: For Specificity Blank and Standard are injected into system, there is no any interference of any peak in blank with the retention time of the analytical peaks [12].




Figure 3: Chromatogram for System suitability

2. Linearity:

Preparation of stock solution: Accurately weigh and transfer 100 mg of Elbasvir and 200 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution)




Figure 4: Calibration graph of Elbasvir




Figure 5: Calibration graph of Grazoprevir

Preparation of Level – I (100 ppm of Elbasvir & 200ppm of Grazoprevir): 0.1 ml of above stock solutions has taken in 10ml of volumetric flask, dilute up to the mark with diluent.

Preparation of Level – II (200ppm of Elbasvir & 400ppm of Grazoprevir):

0.2 ml of above stock solutions has taken in 10ml of volumetric flask, dilute up to the mark with diluent.

Preparation of Level – III (300 ppm of Elbasvir & 600ppm of Grazoprevir):

0.3 ml of above stock solutions has taken in 10ml of volumetric flask, dilute up to the mark with diluent.

Preparation of Level – IV (400ppm of Elbasvir & 800ppm of Grazoprevir):

0.4 ml of above stock solutions has taken in 10ml of volumetric flask, dilute up to the mark with diluent

Preparation of Level – V (500ppm of Elbasvir & 1000ppm of Grazoprevir):

0.5 ml of above stock solutions has taken in 10ml of volumetric flask, dilute up to the mark with diluents[13]

Procedure:

Inject each level into the chromatographic system and measure the peak area. Plot a graph of peak area versus concentration (on X-axis concentration and on Y-axis Peak area) and calculate the correlation coefficient.

3. Precision

International Journal of Chemistry and Pharmaceutical Sciences

Preparation of stock solution:

Accurately weigh and transfer 100 mg of Elbasvir and 200 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluent.

Procedure:

The standard solution was injected for six times and measured the area for all six. Injections in HPLC. The %RSD for the area of six replicate injections was found to be within the specified limits [14].

4. Intermediate Precision/Ruggedness

To evaluate the intermediate precision (also known as Ruggedness) of the method, Precision was performed on different day.

Preparation of stock solution:

Accurately weigh and transfer 100 mg of Elbasvir and 200 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluents [15].

Procedure:

The standard solutions prepared in the precision was injected on the other day, for six times and measured the area for all six injections in HPLC. The %RSD for the area of six replicate injections was found to be within the specified limits.

5. Accuracy

Preparation of Standard stock solution:

Accurately weigh and transfer 100 mg of Elbasvir and 200 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluents [16].

Preparation Sample solutions:

For preparation of 50% solution (With respect to target Assay concentration): Accurately weigh and transfer 50 mg of Elbasvir and 100 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluents [17].

For preparation of 100% solution (With respect to target Assay concentration):

Accurately weigh and transfer 100 mg of Elbasvir and 200 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluent.

For preparation of 150% solution (With respect to target Assay concentration):

Accurately weigh and transfer 150 mg of Elbasvir and 300 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluents [18].

Procedure:

Inject the standard solution, Accuracy -50%, Accuracy -100% and Accuracy -150% solutions.

Calculate the Amount found and Amount added for Elbasvir & Grazoprevir and calculate the individual recovery and mean recovery values.

6. Limit of Detection: (for Elbasvir)

Preparation of 300µg/ml solution:

Accurately weigh and transfer 100 mg of Elbasvir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution)

Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluents [19].

Preparation of 0.99 µg/ml solution:

Further pipette 0.33ml of the above stock solution into a 10ml volumetric flask and dilute up to the mark with diluent.

Limit of Detection: (for Grazoprevir)

Preparation of 600 µg/ml solution:

Accurately weigh and transfer 200 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluents [20].

Preparation of 12.0 µg/ml solution:

Further pipette 0.2ml of the above stock solution into a 10ml volumetric flask and dilute up to the mark with diluent

7. Limit of Quantification

Preparation of 300 µg/ml solution:

Accurately weigh and transfer 100 mg of Elbasvir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluent.

Preparation of 33.0 µg/ml solution: Further pipette 1.1 ml of the above stock solution into a 10ml volumetric flask and dilute up to the mark with diluents.

Preparation of 600µg/ml solution:

Accurately weigh and transfer 200 mg of Grazoprevir working standard into a 10 ml clean dry volumetric flask add about 7 mL of Diluent and sonicate to dissolve it completely and make volume up to the mark with the same solvent. (Stock solution). Further pipette 0.3 ml of the above stock solutions into a 10ml volumetric flask and dilute up to the mark with diluent.

Preparation of 37.8 µg/ml solution:

Further pipette 0.63ml of the above stock solution into a 10ml volumetric flask and dilute up to the mark with diluent.

8. Robustness

As part of the Robustness, deliberate change in the Flow rate, Mobile Phase composition, Temperature Variation was made to evaluate the impact on the method.

A. The flow rate was varied at 0.9 ml/min to 1.1ml/min.

Standard solution 300 ppm of Elbasvir & 600ppm of Grazoprevir was prepared and analysed using the varied flow rates along with method flow rate. On evaluation of the above results, it can be concluded that the variation in flow rate affected the method significantly. Hence it indicates that the method is robust even by change in the flow rate $\pm 10\%$.

B. The Organic composition in the Mobile phase was varied from $\pm 10\%$.

Standard solution 300ppm of Elbasvir & 600 ppm of Grazoprevir was prepared and analysed using the varied Mobile phase composition along with the actual mobile phase composition in the method. On evaluation of the above results, it can be concluded that the variation in 10%. Organic composition in the mobile phase affected the method significantly. Hence it

Indicates that the method is robust even by change in the Mobile phase ± 10

9. Degradation Studies

The International Conference on Harmonization (ICH) guideline entitled stability testing of new drug substances and products requires that stress testing be carried out to elucidate the inherent stability characteristics of the active substance. The aim of this work was to perform the stress degradation studies on the Elbasvir and Grazoprevir using the proposed method.

Preparation of stock:

Accurately weigh 10 tablets crush in mortar and pestle and transfer equivalent to 100 mg of Elbasvir and 200 mg Grazoprevir in sample into a 10mL clean dry volumetric flask add about 7mL of Diluent and sonicate it up to 30 mins to dissolve it completely and make volume up to the mark with the same solvent. Then it is Filtered through 0.44 micron Injection filter. (Stock solution).

Hydrolytic degradation under acidic condition

Pipette 0.3 ml of above solution into a 10ml volumetric flask and 3 ml of 0.1N HCl was added. Then, the volumetric flask was kept at 60°C for 24 hours and then neutralized with 0.1 N NaOH and make up to 10ml with diluent. Filter the solution with 0.44 microns syringe filters and place in vials.

Hydrolytic degradation under alkaline condition

Pipette 0.3 ml of above solution into a 10ml volumetric and add 3ml of 0.1N NaOH was added in 10ml of volumetric flask. Then, the volumetric flask was kept at 60°C for 24 hours and then neutralized with 0.1N HCl and make up to 10ml with diluent. Filter the solution with 0.44 microns syringe filters and place in vials [21].

Thermal induced degradation

Elbasvir and Grazoprevire sample was taken in petridish and kept in Hot air oven at 110⁰ C fo 3 hours. Then the

sample was taken and diluted with diluents and injected into HPLC and analyzed.

Oxidative degradation

Pipette 0.3 ml above stock solution into a 10ml volumetric flask and 1ml of 30% w/v of hydrogen peroxide added in 10 ml of volumetric flask and the volume was made up to the mark with diluent. The volumetric flask was then kept

at room temperature for 15 min. Filter the solution with 0.45 microns syringe filters and place in vials.

Photo degradation:

Pipette 0.3 ml above stock solution into a 10ml volumetric flask and expose to sunlight for 24hrs and the volume was made up to the mark with diluent. Filter the solution with 0.45 microns syringe filters and place in vials

Table 1: Area of different concentration of Elbasvir

S. No	Linearity Level	Concentration	Area
1	I	100	16472
2	II	200	32577
3	III	300	47931
4	IV	400	61145
5	V	500	76795
Correlation Coefficient			0.999

Table 2: Area of different concentration of Grazoprevir

S. No	Linearity Level	Concentration	Area
1	I	200	32441
2	II	400	67728
3	III	600	100630
4	IV	800	134448
5	V	1000	172463
Correlation Coefficient			0.999

Table 3: Accuracy (recovery) data for Elbasvir

%Concentration (at specification Level)	Area	Amount Added (mg)	Amount Found (mg)	% Recovery	Mean Recovery
50%	24492	50	49.95	99.89	100.11
100%	48488	100	99.70	99.70	
150%	73486	150	151.10	100.73	

Table 4: Accuracy (recovery) data for Grazoprevir

%Concentration (at specification Level)	Area	Amount Added (mg)	Amount Found (mg)	% Recovery	Mean Recovery
50%	49702	100	100.52	100.52	100.38
100%	99405	200	201.04	100.52	
150%	148499	300	300.34	100.11	

Table 5: Results of Precision for Elbasvir and Grazoprevir

Injection	Area for Elbasvir	Area for Grazoprevir
Injection-1	48997	98709
Injection-2	48348	98962
Injection-3	48957	98700
Injection-4	48487	98687
Injection-5	48674	98901
Injection-6	48691	98960
Average	48692.3	98819.8
Standard Deviation	254.5	134.7
%RSD	0.5	0.1

Table 6: Results of Intermediate precision for Elbasvir and Grazoprevir

Injection	Area for Elbasvir	Area for Grazoprevir
Injection-1	48673	98783
Injection-2	48720	98674
Injection-3	48793	98647
Injection-4	48657	98359

Injection-5	48082	98747
Injection-6	48956	98911
Average	48646.8	98686.8
Standard Deviation	297.4	185.7
%RSD	0.6	0.2

Table 7: Results for variation in flow for Grazoprevir

S. No	Flow Rate (ml/min)	System Suitability Results	
		USP Plate Count	USP Tailing
1	0.9	3013.80	1.4
2	1.0	2935.56	1.42
3	1.1	2845.18	1.43

Table 8: Results for variation in flow for Grazoprevir

S. No	Flow Rate (ml/min)	System Suitability Results		
		USP Plate Count	USP Tailing	USP Resolution
1	0.9	4951.17	1.46	6.64
2	1.0	4800.53	1.46	6.50
3	1.1	4596.34	1.42	6.34

Table 9: System suitability results for Elbasvir

S. No	Change in Organic Composition in the Mobile Phase	System Suitability Results	
		USP Plate Count	USP Tailing
1	10% less	3013	1.1
2	*Actual	2935.56	1.42
3	10% more	2841.98	1.44

Table 10: System suitability results for Grazoprevir

S. No	Change in Organic Composition in the Mobile Phase	System Suitability Results		
		USP Plate Count	USP Tailing	USP Resolution
1	10% less	4751	1.45	6.64
2	*Actual	4800.53	1.46	6.50
3	10% more	4160.01	1.42	4.04

Table 11: Limit of Detection for Elbasvir and Grazoprevir

Drug name	Baseline noise(μ V)	Signal obtained (μ V)	S/N ratio
Elbasvir	58	172	2.97
Grazoprevir	558	179	3.09

Table 12: Limit of Quantification for Elbasvir and Grazoprevir

Drug name	Baseline noise(μ V)	Signal obtained (μ V)	S/N ratio
Elbasvir	58	575	9.91
Grazoprevir	58	584	10.07

Table 13: Results of degradation parameters

Sample Name	Elbasvir	
	Area	% Degraded
Standard	48538	
Acid	44953	7.39
Base	47130	2.90
Peroxide	45836	5.57
Thermal	44789	7.72
Photo	46184	4.85
Sample Name	Grazoprevir	
	Area	% Degraded
Standard	98691	
Acid	95659	3.07

Base	92821	5.95
Peroxide	97498	1.21
Thermal	93151	5.61
Photo	96389	2.33

4. Conclusion

The estimation of Elbosvir and Grazopravir was done by RP-HPLC. The assay of Elbosvir and Grazopravir was performed with tablets and the % assay was found to be 100.18 and 100.00 which shows that the method is useful for routine analysis. The linearity of Elbosvir and Grazopravir was found to be linear with a correlation coefficient of 0.999 and 0.999, which shows that the method is capable of producing good sensitivity. The acceptance criteria of precision is RSD should be not more than 2.0% and the method show precision 0.5 and 0.1 for Elbosvir and Grazopravir which shows that the method is precise. The acceptance criteria of intermediate precision is RSD should be not more than 2.0% and the method show precision 0.6 and 0.2 for Elbosvir and Grazopravir which shows that the method is repeatable when performed in different days also. The accuracy limit is the percentage recovery should be in the range of 97.0% - 103.0%. The total recovery was found to be 100.11% and 100.38% for Elbosvir and Grazopravir. The validation of developed method shows that the accuracy is well within the limit, which shows that the method is capable of showing good accuracy and reproducibility. The acceptance criteria for LOD and LOQ is 3 and 10. The LOD and LOQ for Elbosvir was found to be 2.97 and 9.91 and LOD and LOQ for Grazopravir was found to be 3.09 and 10.07. The robustness limit for mobile phase variation and flow rate variation are well within the limit, which shows that the method is having good system suitability and precision under given set of conditions.

5. References

- [1] Becket and Stenlake, Practical pharmaceutical chemistry, part 24th edition CBS publications and distributors, 2005.
- [2] P.D. Sethi, HPLC quantitative analysis of pharmaceutical formulations CBS publications and distributors, 1st edition, 2001.
- [3] B.K Sharma, instrumental method of chemical analysis, 23rd edition, goal publishers 2004.
- [4] Validating chromatographic methods, David M. Bliesner. International conference on harmonization: ICH Q 2 (R1) Validation of Analytical Procedures: Text and Methodology 1995.
- [5] Skoog D A, West D M, Holler FJ: Introduction of analytical chemistry. Sounder college of publishing, Harcourt Brace college publishers. (1994), PP 1-5.
- [6] ZEPATIER™ (elbasvir and grazoprevir) Tablets, for Oral Use. Full Prescribing Information" (PDF). Merck Sharp & Dohme Corp., a subsidiary of Merck & Co., Inc. Retrieved 31 January 2016.
- [7] FDA approves Zepatier for treatment of chronic hepatitis C genotypes 1 and 4" U.S. Food and Drug Administration. Retrieved 31 January 2016.
- [8] Lindsay Sandy, HPLC by open learning, (1991) 30-45.
- [9] Lough WJ, Wainer IWW. HPLC fundamental principles and practices, (1991) 52- 67.
- [10] Krstulovic AM, Brown PR. Reversed-Phase High Performance Liquid Chromatography: Theory, Practice and Biomedical Applications, Wiley, New York, (1982).
- [11] U.S. FDA, Title 21 of the U.S. Code of Federal Regulations: 21 CFR 211- Current good manufacturing practice for finished pharmaceuticals.
- [12] U.S. FDA - Guidance for Industry (draft) Analytical Procedures and Methods Validation: Chemistry, Manufacturing, and Controls and Documentation, (2000).
- [13] ISO/IEC 17025, General requirements for the competence of testing and calibration laboratories, (2005).
- [14] International Conference on Harmonization (ICH) of Technical Requirements for the Registration of Pharmaceuticals for Human Use, Validation of analytical procedures: definitions and terminology, Q2A, Geneva (1996).
- [15] International Conference on Harmonization (ICH) of Technical Requirements for the Registration of Pharmaceuticals for Human Use, Validation of analytical procedures: Methodology, Q2B, Geneva (1996).
- [16] U.S. EPA, Guidance for methods development and methods validation for the Resource Conservation and Recovery Act (RCRA) Program, Washington, D.C. (1995).
- [17] General Chapter 1225, Validation of compendial methods, United States Pharmacopeia 30, National Formulary 25, Rockville, Md., USA, The United States Pharmacopeial Convention, Inc., (2007).